

AUG MEDICAL Video Laryngoscope

Application:

Difficult airway intubation (CL from I ~ IV grade), Traumatic airway intubation,
Intubation for cervical activity
Limited patients
Replace routine laryngoscope.
Use for anesthesia department, emergency department, ICU,
Cardiology, etc.
Wild rescue, first aid center, ambulance

Clinical Advantages:

Avoid the injury intubation may cause to the laryngeal structure by
getting a optimal view of the laryngeal structure
Unique anti-fog capacity, no need to preheat, save
time in emergency intubation operation
Complete specifications of blades (Pediatric Size/Adult
Size/Difficult Airway Intubation Size optional)
Ultra-big capacity of the battery, can be continually used for 200
minutes.



Ref: VLR-100

